

Overview of the Dolphin Neurostim Professional

The Dolphin Neurostim Professional is specifically designed to relieve pain. Dolphin owes its exceptional outcomes to the extremely accurate point location and the unique, gentle stimulation the unit produces.

Dolphin produces thousands of times less stimulation than traditional TENS with pads/gels. The safety factor and powerful outcomes of this form of stimulation have been thoroughly tested over 30 years by 50,000 health professionals around the world.

Dolphin's stimulation is so gentle, safe, and easy to apply, that any person may now enjoy on-the-spot pain relief at home. Dolphin stimulation is also very powerful, so professionals needing adjunctive clinical pain relief now resort to the Dolphin!

For therapists, Dolphin is also a particularly attractive solution for patients with an aversion to needles, or those with a fear of Hepatitis or AIDS transmission through skin penetration.

About the Dolphin Neurostim Professional Device

The Dolphin is intended for treatment of pain where no medical or surgical intervention is appropriate. The Dolphin is classified as TENS device, and therefore provides a symptomatic treatment of nonproven curative value and as such suppresses the sensation of pain, which may otherwise serve as a protective mechanism. Treatment should be prescribed only after a medical diagnosis of the cause of the pain has been established.

Technical Description of Dolphin Neurostim Professional

Intended use: Transcutaneous Electrical Nerve Stimulator for chronic and acute pain management.

Point location: Metal Tip is moved around the skin until audible tone rises in pitch and the indicator light flickers off, signifying accurate point location.

Stimulation: Pressing the stimulation button activates the stimulation mode and applies electro-stimulation.